

PI 30 Mobile X-Ray

Looking at a new world



GM General
Medical
Italia

GENERAL MEDICAL ITALIA Ltd.

10 Philpot Lane, London, England, EC3M 8AA
Tel. +44 (0) 207 283 7489 - Fax +44 (0) 203 283 4001

www.gmitalia.co.uk - info@gmitalia.co.uk

PI 30



The new X-RAY generation Analogue & Digital Radiography Solution

*Possibility to upgrade from
analogue to digital*

PI 30 is a mobile system used for diagnosis and X-ray examinations in radiography.

It allows to perform X-rays on CR or film by setting the most suitable radiological data according to the interested anatomic area.

The high handiness allows an easy positioning of the unit close to any patient bed with precise movements thanks to the rotation of the column: $\pm 90^\circ$.



**Microprocessor managed control panel
with user-friendly interface**

ELECTRICAL & RADIOLOGICAL FEATURES

Single phase voltage	115 / 230 Vac \pm 10%
Frequency	50/60 Hz
Maximum Power	32 kW
kV range	40-125 (1 kV step)
mA range	50-400 according to kV selection
mAs range	0.1-220 (12,5% steps)
Time range (according to mAs selection)	0,001 – 2,2 s
Connection to standard mains outlet	16 A

GB *New design.*

MAIN FEATURES:

- Maximum power: 32 kW
- Storage of 36 exams (APR)
- Radiographic technique at 2 points
- Cassette Holder (format 35x43): for 5 cassettes
- Remote control device as optional