

**TAU -30**

Mobile C Arm (32kW)

Looking at a new world



**GMI** General Medical Italia

**GENERAL MEDICAL ITALIA Ltd.**

10 Philpot Lane, London, England, EC3M 8AA  
Tel. +44 (0) 207 283 7489 - Fax +44 (0) 203 283 4001  
[www.gmitalia.co.uk](http://www.gmitalia.co.uk) - [info@gmitalia.co.uk](mailto:info@gmitalia.co.uk)

# TAU - 30




The new X-RAY generation  
Digital Radiography Solution

## TAU - 30

**MOBILE 'C'-arm :**  
with FPD 30x30 or 21x21 cm  
equipped with Digital Angiographic  
Memory suitable for all the applications  
having high diagnostic content.

**USER'S INTERFACE:**  
Easy and practical 10" Touch Screen  
control console, installed directly on the  
unit, for all the operative parameters and  
messages of possible faults.



**DISPLAY STATION:**  
The monitor trolley is equipped  
with a 27" High Resolution main monitor  
for work and memory images displaying  
and a dedicated 21,5" Touch Screen  
monitor for the operator.

- HIGHLIGHTS:**
- 30 kW H.V. Generator
  - Digital Memories: 1kx1k,
  - DUAL COOLING System for immediate and effective heat removal,
  - DUAL POWER System: Power reserve system
  - E-MOTION System: C-arm movements, fully motorized.